



Snyllofit™

CONTOURING MASTERX

MAKE
YOUR
BODY
SLIM

CERTIFIED BY



SnyllofitTM SculptX is the most advanced and intensive robotic magnetic muscle stimulator(RMS).

Focused magnetic stimulation passes through all skin and fat in order to effectively stimulate the muscle, providing the most intensive continuous contractions which is ideal for muscle growth, and also inducing apoptosis.

SnyllofitTM SculptX is an electro-magnetic muscle stimulator device designed for aesthetic purpose, having four (4) applicators with higher intensity.

It's cutting-edge technology in non-invasive body contouring, as it not only burns fat, but also builds muscle.

This is key to a toned and healthy physique.

New Technology in Non-invasive Body Slimming



CERTIFIED BY



SnyllofitTM is a 93/42/EEC MDD certified, TGA approved medical device. Thanks to its patented air cooling system, magnetic output intensity of 8 Tesla it is the most powerful system in its class today unrivalled by the competitors.


Four applicators are placed on the target muscle area, such as abs, thighs or buttocks. The applicators then generate intense magnetic stimulations that cause involuntary muscle contractions.

Four Applicators

These contractions trigger the release of free fatty acids, which break down fat deposits and increase muscle tone and strength.


This is very similar to the way of doing a workout.

As working out, the stimulation causes to rebuild and repair muscle tissue, resulting in stronger, firmer muscles.




Non-invasive medical technology that is used for strengthening and re-education of muscles via interaction of the magnetic field with the tissue of the patient.

Contrary to voluntary muscle contractions, the supramaximal contractions are independent of brain function. The magnetic stimulation uses a specific range of frequencies that does not allow muscle relaxation between two consecutive stimulations.



Proven Technology

The magnetic stimulation uses a specific range of frequencies that does not allow muscle relaxation between two consecutive stimulations. The muscle is forced to remain in a contracted state for multiple seconds. When repeatedly exposed to these high load conditions the muscle tissue is stressed and is forced to adapt.



Mechanism of muscle contraction

Recent studies reported that on average 15% - 16% increase in abdominal muscle thickness was observed in treated patients one to two months after magnetic muscle stimulation treatments.




HIEMT

High Intensity Electromagnetic Muscle Trainer

is the medical technology that is used in aesthetic medicine. It uses focused electromagnetic field with a safe intensity level.

Electromagnetic field passes non-invasively through the body and interacts with motor neurons which subsequently trigger supramaximal muscle contractions.


"HIFEM+RF is used as a non-invasive body contouring treatment for the abdomen, thighs and buttocks."



Effect on fat, several recent studies using CT, MRI and ultrasound evaluations have reported approximately 35% reduction in subcutaneous fat layer in patients treated by a magnetic stimulation based device on their abdomen. The principle of cell apoptosis induced by increased concentrations of free fatty acids has been previously observed and demonstrated in numerous research studies.

Session Information

Snyllofit™ SculptX treatment session requires no anaesthesia, incisions or discomfort. In fact, patients are able to sit back and relax, while the device performs the equivalent of up to 50,000 painless crunches or squats in 30 minutes.



Hercules Course

Building Muscle

Treatment	1	2	3	4	5	6
Program	Hypertrophy	Combo2	Hypertrophy	Combo2	Hypertrophy	Combo2

David Course

Burning Fat

Treatment	1	2	3	4	5	6
Program	HIIT	Combo1	Strength	Strength	HIIT	Combo1

Single treatment session takes 30 minutes only, and there is no downtime. 2-3 treatments per week, would be enough for perfect result for most people. Generally 4-6 treatments recommended.



Workout Program

Snyllofit™ SculptX unique workout program consists of a series of an actual physical training program, which has been ideally designed for the maximum stimulation leading to the best result.

Venus Course Building Muscle


Treatment	1	2	3	4	5	6
Program	Hypertrophy	HIIT	Combo2	Hypertrophy	HIIT	Combo2

Aphrodite Course Burning Fat

Treatment	1	2	3	4	5	6
Program	HIIT	Combo 1	Strength	HIIT	Combo2	HIIT

Additionally, **Snyllofit™ SculptX** is an extremely efficient treatment that takes place over the course of 10-14 days, rather than months, and it also produces more toned abs and uptight buttocks within 2-4 weeks.

If you need more information about workout program, please visit www.snyllo.com and send us your enquiry.





Snyllofit™ Advantages



Output Intensity 8-13 Tesla

HIFEM+RF base 8-13 Tesla output enables that magnetic energy could embrace big skeletal muscles of the human body. That means **Snyllofit™** SculptX muscle re-education process is deep enough like normal muscle contractions. Superficial electrical muscle stimulation is incomparable.



Air Cooling Technology

Our optimised air cooling system makes **Snyllofit™** SculptX effective to operate for a very long time without any downtime.

Even further, there is no worry for leaking oils anymore.



Manual Training

Our cutting-edge technology makes it possible to adjust magnetic pulse width through to a very wide range. So you don't have to use set presets, you can design special manual treatment programs for your patients' needs.



No More Extra Supplies, **Snyllofit™** No Consumables



Perfect Air Cooling System makes **Snyllofit™** Even More Stronger



Great Durability - 25,000 Treatment Assurance



2 Year Full Coverage Warranty



Convenient Operation - Easy & Effective Preset / User Interface



Maintenance

Before & After

Even if you aren't in optimal shape, you can experience benefits from **Snyllofit™** for years to come with regular exercise and balanced eating.




M1_Before_Front

M1_After_Front




F1_Before_Front

F1_After_Front



M1_Before_Side

M1_After_Side



M2_Before_Front

M2_After_Front

FAQs

Q A candidate?

A Anybody

Q A downtime?

A Not at all

Q How long?

A 30 minutes
per session

Q How often?

A 2~3 times
a week

Q How many time?

A 4~6 treatment
recommended

Q How long can
I use it?

A No limit

Q How long will
it last?

A Extremely long

Q Additional
supplies?

A No
consumables

Q Warranty?

A 1 year

Q Manual program?

A Easy & Effective
manual mode

- 1 LCD 2 Emergency stop button 3 Reset Switch 4 Applicator (2EA)

Specification

Power	Power Input	240 V 50/60 Hz
	Maximum Power Consumption	2.3 kVA
Magnetic Stimulation	Magnetic induction amplitude(Intensity)	8 Tesla
	Shape of stimulation pulse	Biphasic Wave
	Pulse duration	320 μ s
Dimension	Size (height x width x depth)	570x376x900 (mm)
	Weight	80 kg
	Safety type	B
Country of Origin	Japanese-German Technology	
Manufacturer	SNYLLO INNOVATIONS PVT LTD	





Snyllofit™

Find out more about **Snyllofit™**



🌐 www.snyllo.com

